



s&s simple&smart

TREATMENT CHAIRS

MADE IN ITALY

TREATMENT CHAIRS



TREATMENT CHAIRS

Elegant, functional, with a unique Italian design...

Find out which chair best suits your demanding taste and expertise.

Built to meet the challenge of your expectations.

LINDA NEXT

AESTHETIC MEDICINE

AUDIOLOGY

DERMATOLOGY

DIABETOLOGY

DIGITAL DENTISTRY

ENDOCRINOLOGY

ENT

GENERAL MEDICAL EXAMINATIONS

LASER TREATMENTS

MAXILLO FACIAL

SURGERY

TRICHOLOGY

Strong and flexible, ergonomic and functional, with a lifting capacity up to 300 kg. Wide range of accessories and options to keep up with the most modern and ambitious professionals needs.

Standard features

- Complete symmetry for right – or left-handed professionals
- Electromechanical lifting, lowering and backrest tilting
- Backrest movement synchronized with leg section
- Automatic extensible PU foot section
- Independent Trendelenburg movement
- Automatic built in programs for patient positioning: table, Trendelenburg, rinse, last and zero position
- 3 operator profiles with 3 programmable positions each
- Safety security STOPS on backrest, leg-rest and controls
- SLOW MODE feature for precise movements
- Maximum lifting capacity 300 Kg
- 24 DC linear actuators
- Soft start & soft stop movement
- Joystick control
- VALE headrest (double articulated)
- Standard molded upholstery
- RAL 9003 painting
- Different PU and fabric colors



TECHNICAL DATA

| | | |
|---|--------|--------|
| Minimum height | mm | 500 |
| Maximum height | mm | 875 |
| Backrest height adjustment | mm | 60 |
| Backrest tilting angle | ° | 0°-90° |
| Width | mm | 630 |
| Length | mm | 2000 |
| Net weight approx. | Kg | 140 |
| Maximum lifting capability (patient + unit) | Kg | 300 |
| Lifting raising speed | mm/sec | 22 |
| Trendelenburg | ° | 21,5° |
| Rotation angle | ° | 325° |
| Classification | | I |
| Class | | B |





Options

- Independent leg section
- Swivelling left and right armrests
- Different headrests available according to specific treatment needs
- Memory or handmade soft foam upholsteries (also with handmade contrast stitching)
- Backrest motorized height adjustment
- Wheels
- Rotation system around the vertical axis under the chair base (320°)
- Backup Battery
- Different chair painting colors
- IR remote control
- Bracket for lamp and/or tray connection
- Children pillow
- Leg section foot protection
- Stabilizing plate



Swivelling armrests



Sliding backrest



Sliding backrest



Soft Foam upholstery



IR Remote control



Children pillow

LINDA NEXT

SPECIFIC CONFIGURATIONS

LINDA NEXT chair, thanks to its versatility, is available in several different possible configurations, each conceived according to specific treatment purposes:

- LINDA NEXT AES: for medical aesthetic procedures
 - LINDA NEXT ENT: for otorhinolaryngology and audiology procedures
 - LINDA NEXT SUR: for daily outpatient surgery
-

LINDA NEXT AES

LINDA NEXT AES is the chair configuration setting designed for aesthetic medicine, which can be equipped with the following specific accessories:

- COMFORT wide multi-articulated headrest
- Bracket for lamp and/or tray connection
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Paper roll holder (roll width 45 cm)
- Polyurethane gel handrest
- High-tech EVA CAM medical lamp with integrated 4K camera, mounted on the chair or ceiling
- Articulated arm with large instrument tray





Standard Options

- Independent leg section
- Different headrests available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Backrest motorized height adjustment
- Wheels
- Rotation system around the vertical axis under the chair base (320°)
- Backup Battery
- Different chair painting colors
- IR remote control
- Leg section foot protection
- Stabilizing plate

AES Options

- Bracket for lamp and/or tray connection
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Paper roll holder (roll width 45 cm)
- Polyurethane gel handrest
- High-tech EVA CAM medical lamp with integrated 4K camera, mounted on the chair or ceiling
- Articulated arm with large instrument tray
- PRIMA FLEX compact directional lamp
- Compact instrument tray with bracket
- COMFORT wide multi-articulated headrest
- Göpel leg holders
- Small side instrument tray with chair connection
- Tablet or articulating mirror holder



Polyurethane gel handrest



Surgical rotating armrest



Göpel leg holders



COMFORT headrest



COMFORT headrest



Compact directional lamp and tray

LINDA NEXT ENT

LINDA NEXT ENT is the chair configuration setting designed for otorhinolaryngology procedures and audiology purposes, which can be equipped with the following specific accessories:

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Swivelling left and right armrests
- LIGHT multi-articulated headrest





Standard Options

- Independent leg section
- Different headrests available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Backrest motorized height adjustment
- Backup Battery
- Different chair painting colors
- Bracket for lamp and/or tray connection
- Children pillow
- Leg section foot protection
- Stabilizing plate

ENT Options

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- LIGHT multi-articulated headrest
- Swivelling left and right armrests
- COMFORT wide multi-articulated headrest
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Small side instrument tray with chair connection
- Paper roll holder (roll width 45 cm)



COMFORT headrest



IR Remote control



Base with rotation



LIGHT headrest



Swivelling armrests



Paper roll holder

LINDA NEXT SUR

LINDA NEXT SUR is the chair configuration setting designed for daily outpatient surgery purposes, which includes the following specific accessories:

- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Bracket for lamp and/or tray connection
- Stainless steel adjustable hydraulic serving cart
- LIGHT multi-articulated headrest
- Backup Battery





Standard Options

- Independent leg section
- Different headrests available according to specific treatment needs
- Memory or handmade soft foam upholsteries (also with handmade contrast stitching)
- Backrest motorized height adjustment
- Rotation system around the vertical axis under the chair base (320°)
- Different chair painting colors
- IR remote control
- Children pillow
- Leg section foot protection
- Stabilizing plate

SUR Options

- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Articulated arm with large instrument tray
- Stainless steel adjustable hydraulic serving cart
- COMFORT wide multi-articulated headrest
- LIGHT multi-articulated headrest
- Bracket for lamp and/or tray connection
- Backup Battery
- Wheels
- Paper roll holder (roll width 45 cm)
- Small side instrument tray with chair connection



COMFORT headrest



Paper roll holder



Hydraulic serving cart



Surgical rotating armrest



LIGHT headrest



Sliding backrest

LINDA EVO

AESTHETIC MEDICINE

AUDIOLOGY

DERMATOLOGY

DIABETOLOGY

DIGITAL DENTISTRY

ENDOCRINOLOGY

ENT

GENERAL MEDICAL EXAMINATIONS

LASER TREATMENTS

MAXILLO FACIAL

SURGERY

TRICHOLOGY

LINDA EVO is the ultimate multifunctional versatile smart knee-break chair, which embodies elegant yet compact modern design, functional style and ergonomics.

Perfectly suitable to treat all patients, even children, elderly people, and persons with reduced mobility.

Standard features

- Fully electromechanical upward, downward and tilt movements of the backrest
- Backrest movement compensated, up to horizontal position
- Tilting of the backrest synchronized with the elongation of the leg rest
- High enveloping backrest
- Complete symmetry, usable by both right and left-handed professionals
- Safety stops on backrest and lifting arm
- Automatic program for complete return to the patient lowering position ("zero")
- Three user programmable customized working positions
- Automatic rinse program with return to the last working position
- Low voltage controls: 5 V
- Joystick control, on the base
- VALE headrest (double articulated)
- Soft-start and soft-stop movement



TECHNICAL DATA

| | | |
|--|----|--------|
| Max. height | mm | 880 |
| Min. height | mm | 550 |
| Backrest tilting angle | ° | 0 -76 |
| Seat tilting angle | ° | +7/+18 |
| Width | mm | 615 |
| Total length | mm | 1750 |
| Gross weight | Kg | 110 |
| Net weight | Kg | 80 |
| Maximum lifting capability (patient + connections) | Kg | 135+75 |
| Trendelenburg | ° | +18 |
| Rotation angle | ° | +/-160 |
| Classification | | I |
| Class | | B |



Options

- Swivelling left and right armrests
- Different headrest models available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Children pillow



Adjustable headrest



Swivelling armrests




Children pillow



Soft foam



Joystick control



IR Remote control

LINDA EVO

SPECIFIC CONFIGURATIONS

LINDA EVO is a simple yet versatile knee-break chair that can be configured to suit various specific treatment purposes:

- LINDA EVO AES: for medical aesthetic procedures
 - LINDA EVO ENT: for otorhinolaryngology and audiology procedures
 - LINDA EVO SUR: for daily outpatient surgery
-

LINDA EVO AES

LINDA EVO AES is the chair configuration setting designed for aesthetic medicine, which can be equipped with the following specific accessories:

- COMFORT wide multi-articulated headrest
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Polyurethane gel handrest
- PRIMA FLEX compact directional lamp





Standard Options

- Swivelling left and right armrests
- Different headrest models available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Rotation system around the vertical axis under the chair base (320°)
- Different chair painting colors
- IR remote control

AES Options

- COMFORT wide multi-articulated headrest
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Polyurethane gel handrest
- PRIMA FLEX compact directional lamp
- Stainless steel adjustable hydraulic serving cart
- Paper roll holder (roll width 45 cm)
- Tablet or articulating mirror holder



Hydraulic serving cart



Surgical rotating armrest



COMFORT headrest



Polyurethane gel handrest

LINDA EVO ENT

LINDA EVO ENT is the chair configuration setting designed for Otorhinolaryngology procedures and audiology purposes, which can be equipped with the following specific accessories:

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Swivelling left and right armrests
- ELLE Mobile headrest
- ENT footrest





Standard Options

- Swivelling left and right armrests
- Different headrest models available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Different chair painting colors
- Children pillow

ENT Options


- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- ELLE Mobile double-articulated headrest
- COMFORT wide multi-articulated headrest
- Paper roll holder (roll width 45 cm)
- ENT footrest




IR Remote control



ELLE Mobile headrest




COMFORT headrest



Base with rotation



Swivelling armrests



ENT footrest

LINDA EVO SUR

LINDA EVO SUR is the chair configuration setting designed for daily outpatient surgery purposes, which includes the following specific accessories:

- Swivelling left and right armrests
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Stainless steel adjustable hydraulic serving cart
- COMFORT wide multi-articulated headrest
- LIGHT multi-articulated headrest





Standard Options

- Swivelling left and right armrests
- Different headrest models available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Rotation system around the vertical axis under the chair base (320°)
- IR remote control

SUR Options

- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Stainless steel adjustable hydraulic serving cart
- COMFORT wide multi-articulated headrest
- LIGHT multi-articulated headrest
- Paper roll holder (roll width 45 cm)



Surgical rotating armrest



COMFORT headrest



LIGHT headrest



Hydraulic serving cart

www.simplesmart.it



Registered office: Via Perola 16, 24021 Albino, Bg - Italy

Headquarters: Via San Carlo 33/35, 24029 Vertova, Bg - Italy

Ph: +39 (0)35 5295636 | info@simplesmart.it | www.simplesmart.it

