

AMALGAM FILTER SIMPLE ONE (DDS2011) – TECHNICAL DATA

FEATURES

The filter can be installed with all types of suction motors

Connected directly to the suction system pipe

No electronic components

Easy to install and maintenance-free

ISO 11143:2008 certified

The filter includes connectors and fittings to connect it to most standard-sized suction pipes

OTHER INFORMATION

The filter must be installed between the dental unit and the suction motors

If the suction motor has a dry suction system, install the filter between the dental unit and the air/water separator

The filter can be connected only to one dental chair. If there are two or more dental units connected to the same suction motor, install one filter for each dental chair



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it