

COMPRESSOR AIR PLUS ONE – TECHNICAL DATA

Power	1.1 kW / 1.5 HP
Voltage	110 / 220 V, 50 / 60 Hz
Pressure	7 bar
Air get in	200 L/min
Air delivery	150 L/min (at 5 bar)
Tank capacity	50 L
N. of dental units	2 + 1
Dryer	None
Soundproof cover	None
Noise	68 dB (A)
Dimensions	40 x 75 x 75h cm
Weight	47 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it