

COMPRESSOR AIR PLUS TWO WITH DRYER AND SOUNDPROOF COVER – TECHNICAL DATA

Power	1.8 kW / 2.5 HP
Voltage	110 / 220 V, 50 / 60 Hz
Pressure	7 bar
Air get in	300 L/min
Air delivery	225 L/min (at 5 bar)
Tank capacity	50 L
N. of dental units	3 + 1
Dryer	Yes
Soundproof cover	Yes
Noise	64 dB (A)
Dimensions	76 x 50 x 108h cm
Weight	61 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it