

## COMPRESSOR AIR PRO ONE WITH DRYER – TECHNICAL DATA

Power	1.1 kW / 1.5 HP
Voltage	110 / 220 V, 50 / 60 Hz
Pressure	7 bar
Air get in	100 L/min
Air delivery	70 L/min (at 5 bar)
Tank capacity	24 L
N. of dental units	1
Dryer	Yes
Soundproof cover	None
Noise	66 dB (A)
Dimensions	40 x 52 x 65h cm
Weight	35 kg



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)