

COMPRESSOR AIR PRO TWO WITH DRYER AND SOUNDPROOF COVER – TECHNICAL DATA

Power	1.1 kW / 1.5 HP
Voltage	110 / 220 V, 50 / 60 Hz
Pressure	7 bar
Air get in	150 L/min
Air delivery	100 L/min (at 5 bar)
Tank capacity	24 L
N. of dental units	1 + 1
Dryer	Yes
Soundproof cover	Yes
Noise	62 dB (A)
Dimensions	76 x 67 x 78h cm
Weight	41 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it