

CONTRA-ANGLE QUARK CAL 1:1 WITH LIGHT – TECHNICAL DATA

Max speed	40,000 rpm
Spray system	E-type internal
Speed transmission	1:1 (blue ring)
Light	Fiber optic
Pushbutton chuck	CA burs, Ø 2.35 mm
Sterilization	Up to 134 °C, 2 bar
Weight	63 g
Standard Conformity	ISO 14457, ISO 3964



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it