

## CONTRA-ANGLE S&S BLUE 1:1 – TECHNICAL DATA

Max speed	40,000 rpm
Speed transmission	1:1 (blue ring)
Light	Fiber optic
Spray water pressure	2.0 bar
Spray air pressure	2.0 – 3.0 bar
Spray air flow rate	Min. 1.5 NI/min (at 2.0 bar)
Cooling air flow	≤9.5 NI/min
Pushbutton chuck	Ø 2.35 mm with bur limit stop
Sterilization	At least 3 minutes at 134 °C, 2 bar
Insert	Dental burs in accordance with DIN EN ISO 1797 type 1
Can be attached to	All motors with a connector in accordance with DIN EN ISO 3964



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)