

CONTRA-ANGLE S&S RED 1:5 – TECHNICAL DATA

| | |
|----------------------|--|
| Max speed | 40,000 rpm |
| Speed transmission | 1:5 (red ring) |
| Light | Fiber optic |
| Spray water pressure | 2.0 bar |
| Spray air pressure | 2.0 bar |
| Spray air flow rate | Min. 1.5 NI/min (at 2.0 bar) |
| Cooling air flow | ≤9.5 NI/min |
| Pushbutton chuck | Ø 1.6 mm with bur limit stop |
| Sterilization | At least 3 minutes at 134 °C, 2 bar |
| Insert | Dental burs in accordance with DIN EN ISO 1797 type 3 |
| Can be attached to | All motors with a connector in accordance with DIN EN ISO 3964 |



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it