

INTRAORAL CAMERA DR'S CAM V – TECHNICAL DATA

Resolution	Full HD (1920 x 1080)
Sensor	5 MP sensor CMOS / VCM
Autofocus	Yes
Light	LED, adjustable brightness and temperature
Rotatable	Yes, 280° twist head
Operation voltage	5 V
Connectable to a monitor	Yes, with the standalone M-station
Connectable to a PC	Yes, with the dedicated software
USB cable length	3 m
Weight	53 g



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it