

## PNEUMATIC MICROMOTOR VORTICE – TECHNICAL DATA

Rotation speed	4,000 – 22,000 rpm
Maximum torque	2.7 Ncm
Spray system	E-type internal
Light	None
Air pressure	2.4 – 3.0 bar (34.8 – 43.5 PSI)
Max air consumption	63 NI/min
Sterilization	Up to 134 °C, 2 bar
Weight	96 g
Standard conformity	ISO 14457, ISO 3964, ISO 9168



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)