

## SCALER TITANUS E – TECHNICAL DATA

Connection	EMS® compatible
Frequency	27 – 30 KHz
Light	None
Sterilization	Up to 134 °C, 2 bar
Weight	45 g
Standard conformity	ISO 22374, IEC 80601-2-60, CEI EN 60601-1-2



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)