

SUCTION MOTOR SMART POWER DRY – TECHNICAL DATA

Power	0.41 kW
Voltage	230 V, 50 Hz
Ring type	Dry
Motor speed	2800 rpm
Vacuum level	110 mbar
Flow	1250 L/min
N. of dental units	1
Noise	56 dB (A)
Dimensions	25 x 30 x 40h cm
Weight	20 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it