

SUCTION MOTOR SMART POWER ONE WITH SOUNDPROOF COVER – TECHNICAL DATA

Power	0.55 – 0.63 kW
Voltage	230 V, 50 Hz
Ring type	Semi-dry (humid)
Motor speed	2800 – 3400 rpm
Vacuum level	170 – 200 mbar
Flow	1333 L/min
N. of dental units	1
Noise	50 dB (A)
Dimensions	33 x 57 x 60h cm
Weight	20 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it