

## SUCTION MOTOR SMART POWER TWO WITH INVERTER – TECHNICAL DATA

Power	0.63 kW
Voltage	230 V, 50 Hz
Ring type	Semi-dry (humid)
Motor speed	0 – 3980 rpm
Vacuum level	200 – 370 mbar
Flow	Max 2300 L/min
N. of dental units	1 – 3
Noise	59 dB (A)
Dimensions	41 x 51 x 61h cm
Weight	31 kg



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)