

SUCTION MOTOR SMART POWER TWO – TECHNICAL DATA

Power	0.7 – 0.83 kW
Voltage	230 V, 50 Hz
Ring type	Semi-dry (humid)
Motor speed	2800 – 3400 rpm
Vacuum level	220 – 290 mbar
Flow	2300 L/min
N. of dental units	2
Noise	55 dB (A)
Dimensions	34 x 43 x 54h cm
Weight	20 kg



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it