

## TURBINE BRAVIA L – TECHNICAL DATA

|                     |  |
|---------------------|--|
| Max speed           | 365,000 rpm  |
| Connection          | KaVo MULTIflex® compatible                               |
| Air pressure        | 2.8 ± 0.2 bar  |
| Max air consumption | 53 NI/min  |
| Light               | Fiber optic  |
| Spray port          | 3 – jet  |
| Bur FG              | Ø 1.6 mm (Type3) ISO 1797-1, max length 25 mm ISO 6360-1 |
| Head diameter       | Ø 12.6 mm  |
| Maximum power       | 22 W   |
| Lubrication         | By LubriONCE® grease or SO2000 service oil               |
| Sterilization       | Up to 134 °C, 2 bar                                      |
| Thermodisinfection  | Yes  |
| Maximum noise level | 60 dBA   |
| Weight              | 45 g   |
| Standard conformity | EN ISO 14457, CEI EN 60601-1, ISO 9168                   |



**SIMPLE & SMART SRL,**

**Sede Legale:** Via Perola 16, 24021 Albino, BG - ITALY

**Sede Operativa:** Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

[www.simplesmart.it](http://www.simplesmart.it)