

TURBINE S&S – TECHNICAL DATA

Max speed	390,000 rpm
Connection	KaVo MULTIflex® compatible
Air pressure	2.2 - 3.0 bar
Air consumption	30 - 45 NI/min
Light	Fiber optic
Bur FG	Ø 1.6 mm, max length 25 mm
Sterilization	Up to 134 °C
Maximum noise level	60 dBA



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it