

TURBINE THALIA – TECHNICAL DATA

Max speed	365,000 rpm
Connection	3 or 4 ways
Air pressure	2.6 ± 0.2 bar
Max air consumption	55 NI/min
Light	None
Spray port	1 – jet
Bur FG	Ø 1.6 mm (Type3) ISO 1797-1, max length 25 mm ISO 6360-1
Head diameter	Ø 12.6 mm
Maximum power	23 W
Lubrication	SO2000 service oil
Sterilization	Up to 134 °C, 2 bar
Thermodisinfection	Yes
Maximum noise level	66 dBA
Weight	52 g
Standard conformity	EN ISO 14457, CEI EN 60601-1, ISO 9168



SIMPLE & SMART SRL,

Sede Legale: Via Perola 16, 24021 Albino, BG - ITALY

Sede Operativa: Via San Carlo 33/35, 24029 Vertova, BG - ITALY

VAT: IT04077160168. - Ph: +39 (0)35 5295636

www.simplesmart.it